


IGFBP5 Polyclonal Antibody

Catalog No.	E-AB-66147	Reactivity	H,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of HeLa cells using IGFBP5 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human IGFBP5 (NP_000590.1).
GeneID	3488
Swissprot	P24593
Synonyms	IGFBP5,IBP5

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

IGFBP5 (Insulin Like Growth Factor Binding Protein 5) is a Protein Coding gene. Diseases associated with IGFBP5 include Insulin-Like Growth Factor I and Renal Osteodystrophy. Among its related pathways are Myometrial Relaxation and Contraction Pathways and Senescence and Autophagy in Cancer. Gene Ontology (GO) annotations related to this gene include fibronectin binding and insulin-like growth factor I binding. An important paralog of this gene is IGFBP3.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.