


PGRMC2 Polyclonal Antibody

Catalog No.	E-AB-66726	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of NIH/3T3 cells using PGRMC2 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	A synthetic peptide of human PGRMC2
GeneID	10424
Swissprot	O15173
Synonyms	PGRMC2,DG6,PMBP

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

PGRMC2 (Progesterone Receptor Membrane Component 2) is a Protein Coding gene. Diseases associated with PGRMC2 include Arthrogyrosis, Distal, Type 6 and Arthrogyrosis, Distal, Type 7. Gene Ontology (GO) annotations related to this gene include heme binding and steroid binding. An important paralog of this gene is PGRMC1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.