


HBEGF Polyclonal Antibody


Catalog No.	E-AB-66783	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human colon carcinoma using HBEGF Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human placenta using HBEGF Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human HBEGF (NP_001936.1).
GeneID	1839
Swissprot	Q99075
Synonyms	HBEGF,DTR,DTS,DTSF,HEGFL

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

HBEGF (Heparin Binding EGF Like Growth Factor) is a Protein Coding gene. Diseases associated with HBEGF include Diphtheria and Perinatal Necrotizing Enterocolitis. Among its related pathways are Signaling by PTK6 and RET signaling. Gene Ontology (GO) annotations related to this gene include growth factor activity and epidermal growth factor receptor binding. An important paralog of this gene is AREG.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.