


SLIT1 Polyclonal Antibody

Catalog No. E-AB-66817
Storage Store at -20°C. Avoid freeze / thaw cycles.
Applications WB

Reactivity H
Host Rabbit
Isotype IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of 293T cells using SLIT1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen Recombinant fusion protein of human SLIT1 (NP_003052.2).
GeneID 6585
Swissprot O75093
Synonyms SLIT1,MEGF4,SLIL1,SLIT-1,SLIT3

Product Information

Calculated MW 90kDa/167kDa
Observed MW 168kDa
Buffer PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify Affinity purification
Dilution WB 1:500-1:2000

Background

SLIT1 (Slit Guidance Ligand 1) is a Protein Coding gene. Diseases associated with SLIT1 include Hepatic Adenomas, Familial and Diaphragm Disease. Among its related pathways are Development Slit-Robo signaling and Guidance Cues and Growth Cone Motility. Gene Ontology (GO) annotations related to this gene include calcium ion binding and Roundabout binding. An important paralog of this gene is SLIT2.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.