


SUPT6H Polyclonal Antibody

Catalog No.	E-AB-66819	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U-2 OS cells using SUPT6H Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human SUPT6H (NP_003161.2).
GeneID	6830
Swissprot	Q7KZ85
Synonyms	SUPT6H,SPT6,SPT6H,emb-5,SPT6 homolog

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

SUPT6H (SPT6 Homolog, Histone Chaperone And Transcription Elongation Factor) is a Protein Coding gene. Diseases associated with SUPT6H include Immunodeficiency With Hyper-Igm, Type 2 and Granulomatous Disease, Chronic, Autosomal Recessive, Cytochrome B-Positive, Type I. Among its related pathways are Chromatin Regulation / Acetylation and Gene Expression. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and DNA-binding transcription factor activity.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.