


VCPIP1 Polyclonal Antibody

Catalog No.	E-AB-66859	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using VCPIP1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human VCPIP1 (NP_079330.2).
GeneID	80124
Swissprot	Q96JH7
Synonyms	VCPIP1,DUBA3,VCIP135

Product Information

Calculated MW	134kDa
Observed MW	134kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

VCPIP1 (Valosin Containing Protein Interacting Protein 1) is a Protein Coding gene. Diseases associated with VCPIP1 include Pseudosarcomatous Fibromatosis and Charcot-Marie-Tooth Disease, Axonal, Type 2E. Gene Ontology (GO) annotations related to this gene include thiol-dependent ubiquitin-specific protease activity.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.