


HRNR Polyclonal Antibody

Catalog No. E-AB-66877
Storage Store at -20°C. Avoid freeze / thaw cycles.
Applications IF

Reactivity H
Host Rabbit
Isotype IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of HeLa cells using HRNR Polyclonal Antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen Recombinant fusion protein of human HRNR (NP_001009931.1).
GeneID 388697
Swissprot Q86YZ3
Synonyms HRNR,FLG3,S100A16,S100a18,hornerin

Product Information

Buffer PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify Affinity purification
Dilution IF 1:50-1:200

Background

HRNR (Hornerin) is a Protein Coding gene. Diseases associated with HRNR include Epidermolysis Bullosa Simplex With Nail Dystrophy and Ichthyosis Vulgaris. Among its related pathways are Innate Immune System. Gene Ontology (GO) annotations related to this gene include calcium ion binding. An important paralog of this gene is FLG2.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.