


ARFIP2 Polyclonal Antibody


Catalog No.	E-AB-67055	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human lung cancer using ARFIP2 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human breast cancer using ARFIP2 Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human ARFIP2 (NP_036534.1).
GeneID	23647
Swissprot	P53365
Synonyms	ARFIP2,POR1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

ARFIP2 (ADP Ribosylation Factor Interacting Protein 2) is a Protein Coding gene. Diseases associated with ARFIP2 include Gastric Papillary Adenocarcinoma and Signet Ring Cell Adenocarcinoma. Among its related pathways are G-protein signaling RAC1 in cellular process and Development EGFR signaling via small GTPases. Gene Ontology (GO) annotations related to this gene include GTP binding and Rac GTPase binding. An important paralog of this gene is ARFIP1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.