

FCF1 Polyclonal Antibody

Catalog No. E-AB-67063

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Rat
Immunogen	A synthetic peptide of human FCF1.
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB	1:500-1:2000
IF	1:50-1:200

Data




Western blot analysis of extracts of 293T cells using FCF1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:23kDa
Calculated Mw:23kDa



Immunofluorescence analysis of C6 cells using FCF1 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using FCF1 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

FCF1 (FCF1 RRNA-Processing Protein) is a Protein Coding gene. Diseases associated with FCF1 include Progressive External Ophthalmoplegia With Mitochondrial Dna Deletions, Autosomal Dominant 2. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. An important paralog of this gene is UTP23.