


NUMBL Polyclonal Antibody

Catalog No.	E-AB-67415	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of Mouse kidney cells using NUMBL Polyclonal Antibody.

Immunogen Information

Immunogen	Recombinant fusion protein of human NUMBL (NP_004747.1).
GeneID	9253
Swissprot	Q9Y6R0
Synonyms	NUMBL,CAG3A,CTG3a,NBL,NUMB-R,NUMBLIKE,NUMBR,TNRC23

Product Information

Calculated MW	64kDa
Observed MW	65kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

NUMBL (NUMB Like Endocytic Adaptor Protein) is a Protein Coding gene. Diseases associated with NUMBL include Differentiating Neuroblastoma. Among its related pathways are Notch Signaling Pathway (WikiPathways). An important paralog of this gene is NUMB.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.