


STMN4 Polyclonal Antibody

Catalog No.	E-AB-67534	Reactivity	M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using STMN4 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human STMN4 (NP_110422.2).
GeneID	81551
Swissprot	Q9H169
Synonyms	STMN4,RB3

Product Information

Calculated MW	19kDa/22kDa/23kDa/25kDa
Observed MW	33kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:1000

Background

STMN4 (Stathmin 4) is a Protein Coding gene. Diseases associated with STMN4 include Spastic Paraplegia 36, Autosomal Dominant. Gene Ontology (GO) annotations related to this gene include tubulin binding. An important paralog of this gene is STMN3.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.