


FAM117B Polyclonal Antibody

Catalog No.	E-AB-67582	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using FAM117B Polyclonal Antibody.

Immunogen Information

Immunogen	Recombinant fusion protein of human FAM117B (NP_775782.2).
GeneID	150864
Swissprot	Q6P1L5
Synonyms	FAM117B,ALS2CR13

Product Information

Calculated MW	29kDa/61kDa
Observed MW	65kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

FAM117B (Family With Sequence Similarity 117 Member B) is a Protein Coding gene. Diseases associated with FAM117B include Amyotrophic Lateral Sclerosis 2, Juvenile and Lateral Sclerosis. An important paralog of this gene is GLCCI1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.