


EML2 Polyclonal Antibody

Catalog No.	E-AB-67621	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using EML2 Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human EML2 (NP_036287.1).
GeneID	24139
Swissprot	O95834
Synonyms	EML2,ELP70,EMAP-2,EMAP2

Product Information

Calculated MW	70kDa/86kDa/92kDa
Observed MW	71kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

EML2 (EMAP Like 2) is a Protein Coding gene. Diseases associated with EML2 include Sulfonamide Allergy. Gene Ontology (GO) annotations related to this gene include signaling receptor binding and protein C-terminus binding. An important paralog of this gene is EML1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.