


TDRD3 Polyclonal Antibody

Catalog No.	E-AB-67718	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of U2OS cells using TDRD3 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human TDRD3 (NP_110421.1).
GeneID	81550
Swissprot	Q9H7E2
Synonyms	TDRD3

Product Information

Calculated MW	73kDa/83kDa
Observed MW	73kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

TDRD3 (Tudor Domain Containing 3) is a Protein Coding gene. Diseases associated with TDRD3 include Fragile X Syndrome and Chromosome 22Q11.2 Deletion Syndrome, Distal. Among its related pathways are Chromatin Regulation / Acetylation and Translational Control. Gene Ontology (GO) annotations related to this gene include chromatin binding.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.