


PJA2 Polyclonal Antibody

Catalog No.	E-AB-67755	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U2OS cells using PJA2 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human PJA2 (NP_055634.3).
GeneID	9867
Swissprot	O43164
Synonyms	PJA2,Neurodap1,RNF131

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

PJA2 (Praja Ring Finger Ubiquitin Ligase 2) is a Protein Coding gene. Among its related pathways are Neuroscience and Innate Immune System. Gene Ontology (GO) annotations related to this gene include ligase activity and protein kinase A catalytic subunit binding. An important paralog of this gene is PJA1.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.