


LELP1 Polyclonal Antibody

Catalog No.	E-AB-67827	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using LELP1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human LELP1 (NP_001010857.1).
GeneID	149018
Swissprot	Q5T871
Synonyms	LELP1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

LELP1 (Late Cornified Envelope Like Proline Rich 1) is a Protein Coding gene. Diseases associated with LELP1 include Skin Disease. Among its related pathways are Keratinization and Developmental Biology. An important paralog of this gene is LCE5A.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research
Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.