


DPF1 Polyclonal Antibody

Catalog No.	E-AB-67856	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using DPF1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human DPF1 (NP_001128627.1).
GeneID	8193
Swissprot	Q92782
Synonyms	DPF1,BAF45b,NEUD4,neuro-d4

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

DPF1 (Double PHD Fingers 1) is a Protein Coding gene. Diseases associated with DPF1 include Small-Cell Carcinoma Of The Ovary Of Hypercalcemic Type and Spinal Meningioma. Among its related pathways are Respiratory electron transport, ATP synthesis by chemiosmotic coupling, and heat production by uncoupling proteins. and Gastric cancer. Gene Ontology (GO) annotations related to this gene include nucleic acid binding. An important paralog of this gene is DPF2.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.