


ZBTB25 Polyclonal Antibody

Catalog No.	E-AB-67880	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using ZBTB25 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human ZBTB25 (NP_008908.2).
GeneID	7597
Swissprot	P24278
Synonyms	ZBTB25,C14orf51,KUP,ZNF46

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

ZBTB25 (Zinc Finger And BTB Domain Containing 25) is a Protein Coding gene. Diseases associated with ZBTB25 include Interstitial Nephritis, Karyomegalic. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity. An important paralog of this gene is ZBTB49.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.