

Monoclonal Mouse Anti-Human IgE(C3-domain specific)(Biotin conjugated)

Catalog No. E-AB-1073

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Human IgE C3 domain

Purification Method: Affinity purification

Conjugation: Biotin

Concentration: 0.2 mg/mL

Applications: IHC, WB, ELISA

Dilution: Western Blotting: 1:200-4000(Enhanced chemiluminescent detection)

Western Blotting: 1:200-400(DAB detection)

Immunohistochemistry: 1:100-200

Direct ELISA: 1:1000-5000(TMB detection)

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Buffer: PBS with 1% protective protein, 0.04% Proclin 300 and 50% glycerol, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research
Service for life science

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation