


IL1A Polyclonal Antibody


Catalog No.	E-AB-40407	Reactivity	R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Rat Lymphoid nodule using IL1A Polyclonal Antibody at dilution of 1:200



Immunohistochemistry of paraffin-embedded Rat spleen using IL1A Polyclonal Antibody at dilution of 1:200

Immunogen Information

Immunogen	Recombinant Rat Interleukin-1 alpha protein
GeneID	24493
Swissprot	P16598
Synonyms	BAF, FAF, Hematopoietin 1, IL 1 alpha, IL 1A, IL1 ALPHA, IL1, IL1F1

Product Information

Buffer	PBS with 0.05% Proclin300, 50% glycerol, pH7.3.
Purify	Antigen Affinity Purification
Dilution	IHC 1:50-1:200

Background

proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.