


FGFR3 Polyclonal Antibody

Catalog No.	E-AB-40423	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human kidney using FGFR3 Polyclonal Antibody at dilution of 1:100

Immunogen Information

Immunogen	Recombinant Human Fibroblast growth factor receptor 3 protein
GeneID	2261
Swissprot	P22607
Synonyms	ACH, CD333, CD333 antigen, CEK 2, FGFR 3, FGFR3, MFR 3, SAM 3

Product Information

Buffer	PBS with 0.05% Proclin300, 50% glycerol, pH7.3.
Purify	Antigen Affinity Purification
Dilution	IHC 1:50-1:200

Background

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.