

Anti-Ebola virus EBOV(subtype Bundibugyo, strain Uganda 2007) VP24 Monoclonal Antibody

E-AB-V1105

Application	ELISA	Host	Mouse
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Clone No.	08

Important Note Centrifuge before opening to ensure complete recovery of vial contents.

Product Details

Immunogen	Recombinant EBOV (subtype Bundibugyo, strain Uganda 2007) VP24 Protein (His Tag)
Isotype	IgG2b
Host	Mouse
Clone No.	08
Reactivity	Ebola virus
Dilution	ELISA 1:1000-1:2000
Storage Buffer	0.2 µm filtered solution in PBS
Stability & Storage	Ships on ice packs. Store at -20°C
Description	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified Recombinant EBOV (subtype Bundibugyo, strain Uganda 2007) VP24 Protein (His Tag). And the antibody was puri

Antigen Information

Alternate Names	VP24
------------------------	------

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation