

Anti-HIV-1 p24 Protein(group M, subtype B, strain 92418) Monoclonal Antibody

E-AB-V1225

Application	ELISA	Host	Mouse
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Clone No.	9B6C2B7

Important Note Centrifuge before opening to ensure complete recovery of vial contents.

Product Details

Immunogen	Recombinant HIV-1 p24 Protein (group M, subtype B, strain 92418) (His Tag)
Isotype	IgG1
Host	Mouse
Clone No.	9B6C2B7
Reactivity	HIV
Dilution	ELISA 1:1000-1:2000
Storage Buffer	0.2 µm filtered solution in PBS
Stability & Storage	Ships on ice packs. Store at -20°C
Description	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified Recombinant HIV-1 p24 Protein (group M, subtype B, strain 92418) (His Tag). And the antibody was purified b

Antigen Information

Alternate Names	Gag-p24
------------------------	---------

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation