

## Anti-HIV-1 p24 Protein(group M, subtype B, strain 92418) Monoclonal Antibody

E-AB-V1228

<b>Application</b>	ELISA	<b>Host</b>	Rabbit
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.	<b>Clone No.</b>	002

**Important Note** Centrifuge before opening to ensure complete recovery of vial contents.

### Product Details

<b>Immunogen</b>	Recombinant HIV-1 p24 Protein (group M, subtype B, strain 92418) (His Tag)
<b>Isotype</b>	IgG
<b>Host</b>	Rabbit
<b>Clone No.</b>	002
<b>Reactivity</b>	HIV
<b>Dilution</b>	ELISA 1:5000-1:10000
<b>Storage Buffer</b>	0.2 µm filtered solution in PBS
<b>Stability &amp; Storage</b>	Ships on ice packs. Store at -20°C
<b>Description</b>	This antibody was obtained from a rabbit immunized with purified Recombinant HIV-1 p24 Protein (group M, subtype B, strain 92418) (His Tag). And the antibody was purified by Protein A Affinity.

### Antigen Information

<b>Alternate Names</b>	Gag-p24
------------------------	---------

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
Service for life science**

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation