

Anti-Rift Valley fever virus(TAN/Dod-002/07) RVFV-L Polyclonal Antibody

E-AB-V1312

Application	ELISA	Host	Rabbit
Storage	Store at -20°C. Avoid freeze / thaw cycles.		

Important Note Centrifuge before opening to ensure complete recovery of vial contents.

Product Details

Immunogen	Recombinant Rift Valley fever virus (TAN/Dod-002/07) RVFV-L Protein (His Tag)
Isotype	IgG
Host	Rabbit
Reactivity	RVFV
Dilution	ELISA 1:1000-1:2000
Storage Buffer	0.2 µm filtered solution in PBS
Stability & Storage	Ships on ice packs. Store at -20°C
Description	This antibody was produced in rabbits immunized with purified Recombinant Rift Valley fever virus (TAN/Dod-002/07) RVFV-L Protein (His Tag). And the antibody was purified by Protein A affinity chromatography..

Antigen Information

Alternate Names	RVFV-L
------------------------	--------

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation