


SNRPF Polyclonal Antibody

Catalog No.	E-AB-63643	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using SNRPF Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human SNRPF (NP_003086.1).
GeneID	6636
Swissprot	P62306
Synonyms	SNRPF,SMF,Sm-F,snRNP-F

Product Information

Calculated MW	9kDa
Observed MW	12kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

SNRPF (Small Nuclear Ribonucleoprotein Polypeptide F) is a Protein Coding gene. Among its related pathways are Transport of the SLBP independent Mature mRNA and mRNA Splicing - Minor Pathway. Gene Ontology (GO) annotations related to this gene include RNA binding. An important paralog of this gene is LSM6.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.