


REG3A Polyclonal Antibody


Catalog No.	E-AB-64197	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Images



Immunohistochemistry of paraffin-embedded Mouse brain using REG3A Polyclonal Antibody



Immunohistochemistry of paraffin-embedded Rat lung using REG3A Polyclonal Antibody



Immunohistochemistry of paraffin-embedded Rat heart using REG3A Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human REG3A (NP_002571.1).
GeneID	5068
Swissprot	Q06141
Synonyms	REG3A,HIP,HIP/PAP,INGAP,PAP,PAP-H,PAP1,PBCGF,REG-III,REG3

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

This gene encodes a pancreatic secretory protein that may be involved in cell proliferation or differentiation. It has similarity to the C-type lectin superfamily. The enhanced expression of this gene is observed during pancreatic inflammation and liver carcinogenesis. The mature protein also functions as an antimicrobial protein with antibacterial activity. Alternate splicing results in multiple transcript variants that encode the same protein.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.