


PHOSPHO2 Polyclonal Antibody

Catalog No.	E-AB-64336	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using PHOSPHO2 Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human PHOSPHO2 (NP_001008489.1).
GeneID	493911
Swissprot	Q8TCD6
Synonyms	PHOSPHO2,phosphatase,orphan 2

Product Information

Calculated MW	27kDa
Observed MW	27kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:1000

Background

PHOSPHO2 (Phosphatase, Orphan 2) is a Protein Coding gene. Diseases associated with PHOSPHO2 include Robinow Syndrome, Autosomal Dominant 3 and Tylosis With Esophageal Cancer. Among its related pathways are Metabolism and Vitamin B6 metabolism. Gene Ontology (GO) annotations related to this gene include phosphatase activity and pyridoxal phosphatase activity. An important paralog of this gene is PHOSPHO1.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.