


KLF15 Polyclonal Antibody

Catalog No.	E-AB-64625	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using KLF15 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human KLF15 (NP_054798.1).
GeneID	28999
Swissprot	Q9UIH9
Synonyms	KLF15,KKLF

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

KLF15 (Kruppel Like Factor 15) is a Protein Coding gene. Diseases associated with KLF15 include Glomerular Disease and Aortic Aneurysm. Among its related pathways are Adipogenesis and BMAL1-CLOCK,NPAS2 activates circadian gene expression. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity and proximal promoter DNA-binding transcription activator activity, RNA polymerase II-specific. An important paralog of this gene is KLF5.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.