


URM1 Polyclonal Antibody

Catalog No.	E-AB-64834	Reactivity	M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of Rat eye using URM1 Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human URM1 (NP_001252511.1).
GeneID	81605
Swissprot	Q9BTM9
Synonyms	URM1,C9orf74

Product Information

Calculated MW	7kDa/11kDa/15kDa
Observed MW	11kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

URM1 (Ubiquitin Related Modifier 1) is a Protein Coding gene. Diseases associated with URM1 include Benign Epilepsy With Centrotemporal Spikes. Among its related pathways are Gene Expression and tRNA processing. Gene Ontology (GO) annotations related to this gene include sulfurtransferase activity.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.