


ATRNL1 Polyclonal Antibody

Catalog No.	E-AB-65138	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of Mouse liver using ATRNL1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human ATRNL1 (NP_001263211.1).
GeneID	26033
Swissprot	Q5VV63
Synonyms	ATRNL1,ALP,bA338L11.1,bA454H24.1

Product Information

Calculated MW	51kDa/152kDa
Observed MW	153kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

ATRNL1 (Attractin Like 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include carbohydrate binding. An important paralog of this gene is ATRN.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.