


TRNAU1AP Polyclonal Antibody

Catalog No.	E-AB-65327	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using TRNAU1AP Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human TRNAU1AP (NP_060316.1).
GeneID	54952
Swissprot	Q9NX07
Synonyms	TRNAU1AP,PRO1902,SECP43,TRSPAP1

Product Information

Calculated MW	20kDa/32kDa
Observed MW	36kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:200-1:2000

Background

TRNAU1AP (TRNA Selenocysteine 1 Associated Protein 1) is a Protein Coding gene. Diseases associated with TRNAU1AP include Neuronopathy, Distal Hereditary Motor, Type Viib and Combined D-2- And L-2-Hydroxyglutaric Aciduria. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is TIA1.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.