


LARP1 Polyclonal Antibody

Catalog No.	E-AB-65483	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using LARP1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant Protein of human LARP1
GeneID	23367
Swissprot	Q6PKG0
Synonyms	LARP1,LARP

Product Information

Calculated MW	116kDa/123kDa
Observed MW	163kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

LARP1 (La Ribonucleoprotein 1, Translational Regulator) is a Protein Coding gene. Diseases associated with LARP1 include Hepatocellular Carcinoma. Among its related pathways are Translational Control. Gene Ontology (GO) annotations related to this gene include translation initiation factor binding. An important paralog of this gene is LARP1B.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.