


NEDD1 Polyclonal Antibody

Catalog No.	E-AB-62567	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of L929 cells using NEDD1 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human NEDD1 (NP_690869.1).
GeneID	121441
Swissprot	Q8NHV4
Synonyms	NEDD1,GCP-WD,TUBGCP7

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

NEDD1 (NEDD1 Gamma-Tubulin Ring Complex Targeting Factor) is a Protein Coding gene. Diseases associated with NEDD1 include Seckel Syndrome 7 and Orofaciodigital Syndrome Ix. Among its related pathways are Organelle biogenesis and maintenance and Regulation of PLK1 Activity at G2/M Transition. An important paralog of this gene is WDR62.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.