


FZR1 Polyclonal Antibody

Catalog No.	E-AB-61804	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of HeLa cells using FZR1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human FZR1 (NP_057347.2).
GeneID	51343
Swissprot	Q9UM11
Synonyms	FZR1,CDC20C,CDH1,FZR,FZR2,HCDH,HCDH1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

FZR1 (Fizzy And Cell Division Cycle 20 Related 1) is a Protein Coding gene. Diseases associated with FZR1 include Breast Cancer. Among its related pathways are Regulation of activated PAK-2p34 by proteasome mediated degradation and Cell cycle_Role of APC in cell cycle regulation. An important paralog of this gene is CDC20.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.