


NOLC1 Polyclonal Antibody


Catalog No.	E-AB-63828	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using NOLC1 Polyclonal Antibody at dilution of 1:1000.



Immunofluorescence analysis of U-2 OS cells using NOLC1 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human NOLC1 (NP_004732.2).
GeneID	9221
Swissprot	Q14978
Synonyms	NOLC1, NOPP130, NOPP140, NS5ATP13, P130

Product Information

Calculated MW	73kDa/74kDa
Observed MW	120kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:1000-1:2000 IF 1:50-1:200

Background

NOLC1 (Nucleolar And Coiled-Body Phosphoprotein 1) is a Protein Coding gene. Diseases associated with NOLC1 include Lissencephaly and Treacher Collins Syndrome 1. Gene Ontology (GO) annotations related to this gene include GTP binding.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.