


## DTYMK Polyclonal Antibody

<b>Catalog No.</b>	E-AB-61913	<b>Reactivity</b>	H,M,R
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.	<b>Host</b>	Rabbit
<b>Applications</b>	IF	<b>Isotype</b>	IgG

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Images



Immunofluorescence analysis of A549 cells using DTYMK Polyclonal Antibody

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein of human DTYMK (NP_036277.2).
<b>GeneID</b>	1841
<b>Swissprot</b>	P23919
<b>Synonyms</b>	DTYMK,CDC8,PP3731,TMPK,TYMK

### Product Information

<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purify</b>	Affinity purification
<b>Dilution</b>	IF 1:50-1:200

### Background

DTYMK (Deoxythymidylate Kinase) is a Protein Coding gene. Among its related pathways are superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis and Purine metabolism (KEGG). GO annotations related to this gene include kinase activity and thymidylate kinase activity.

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.