


TFPT Polyclonal Antibody

Catalog No.	E-AB-61981	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U2OS cells using TFPT Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human TFPT (NP_037474.1).
GeneID	29844
Swissprot	P0C1Z6
Synonyms	TFPT,FB1,INO80F,amida

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

TFPT (TCF3 Fusion Partner) is a Protein Coding gene. Diseases associated with TFPT include Childhood Leukemia and Retinitis Pigmentosa. Among its related pathways are Metabolism of proteins and Regulation of activated PAK-2p34 by proteasome mediated degradation. Gene Ontology (GO) annotations related to this gene include protein heterodimerization activity and protein kinase binding.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.