


EAF2 Polyclonal Antibody

Catalog No.	E-AB-61538	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A-549 cells using EAF2 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human EAF2 (NP_060926.2).
GeneID	55840
Swissprot	Q96CJ1
Synonyms	EAF2,BM040,TRAITS,U19

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

EAF2 (ELL Associated Factor 2) is a Protein Coding gene. Diseases associated with EAF2 include Eaf. Among its related pathways are Chromatin Regulation / Acetylation and Gene Expression. Gene Ontology (GO) annotations related to this gene include RNA polymerase II regulatory region sequence-specific DNA binding and DNA-binding transcription activator activity, RNA polymerase II-specific. An important paralog of this gene is EAF1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.