


U2AF1L4 Polyclonal Antibody

Catalog No.	E-AB-61557	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A-549 cells using U2AF1L4 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human U2AF1L4 (NP_659424.2).
GeneID	199746
Swissprot	Q8WU68
Synonyms	U2AF1L4,U2AF1-RS3,U2AF1L3,U2AF1L3V1,U2AF1RS3,U2af26

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

U2AF1L4 (U2 Small Nuclear RNA Auxiliary Factor 1 Like 4) is a Protein Coding gene. Diseases associated with U2AF1L4 include Myelophthistic Anemia and Atypical Chronic Myeloid Leukemia. Among its related pathways are Cleavage of Growing Transcript in the Termination Region and Transport of Mature Transcript to Cytoplasm. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is U2AF1.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.