


NGFRAP1 Polyclonal Antibody

Catalog No.	E-AB-61628	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of MCF-7 cells using NGFRAP1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human NGFRAP1 (NP_996798.1).
GeneID	27018
Swissprot	Q00994
Synonyms	BEX3,Bex,DXS6984E,HGR74,NADE,NGFRAP1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

Protein BEX3 is a protein that in humans is encoded by the NGFRAP1 gene. May be a signaling adapter molecule involved in p75^{NTR}-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death. By similarity. May play an important role in the pathogenesis of neurogenetic diseases.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.