


POU6F1 Polyclonal Antibody

Catalog No.	E-AB-61631	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A-549 cells using POU6F1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human POU6F1 (NP_002693.3).
GeneID	5463
Swissprot	Q14863
Synonyms	POU6F1,BRN5,MPOU,TCFB1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

POU6F1 (POU Class 6 Homeobox 1) is a Protein Coding gene. Diseases associated with POU6F1 include Clear Cell Adenocarcinoma Of The Ovary and Clear Cell Adenocarcinoma. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity. An important paralog of this gene is POU6F2.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.