


SMYD3 Polyclonal Antibody

Catalog No. E-AB-61640
Storage Store at -20°C. Avoid freeze / thaw cycles.
Applications IF

Reactivity H,M
Host Rabbit
Isotype IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A-549 cells using SMYD3 Polyclonal Antibody

Immunogen Information

Immunogen Recombinant fusion protein of human SMYD3 (NP_073580.1).
GeneID 64754
Swissprot Q9H7B4
Synonyms SMYD3,KMT3E,ZMYND1,ZNFN3A1,bA74P14.1

Product Information

Buffer PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify Affinity purification
Dilution IF 1:50-1:100

Background

This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.