


SPIN2B Polyclonal Antibody

Catalog No.	E-AB-62602	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human kidney using SPIN2B Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human SPIN2B (NP_001006682.1).
GeneID	474343
Swissprot	Q9BPZ2
Synonyms	SPIN2B,SPIN-2,SPIN-2B,SPIN2_duplicate,TDRD26,dJ323P24.2

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

SPIN2B (Spindlin Family Member 2B) is a Protein Coding gene. Diseases associated with SPIN2B include Rapp-Hodgkin Syndrome and Myasthenic Syndrome, Congenital, 13. An important paralog of this gene is SPIN2A.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.