


OSGEP Polyclonal Antibody

Catalog No.	E-AB-61731	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using OSGEP Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human OSGEP (NP_060277.1).
GeneID	55644
Swissprot	Q9NPF4
Synonyms	OSGEP,GCPL1,KAE1,OSGEP1,PRSMG1,TCS3

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

OSGEP (O-Sialoglycoprotein Endopeptidase) is a Protein Coding gene. Diseases associated with OSGEP include Galloway-Mowat Syndrome 3 and Galloway-Mowat Syndrome. Among its related pathways are Gene Expression and tRNA processing. Gene Ontology (GO) annotations related to this gene include transferase activity, transferring acyl groups other than amino-acyl groups and transferase activity, transferring acyl groups. An important paralog of this gene is OSGEPL1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.