


PPP3R2 Polyclonal Antibody

| | | | |
|---------------------|---|-------------------|--------|
| Catalog No. | E-AB-61748 | Reactivity | H,M,R |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. | Host | Rabbit |
| Applications | WB | Isotype | IgG |

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using PPP3R2 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

| | |
|------------------|---|
| Immunogen | Recombinant fusion protein of human PPP3R2 (NP_671709.1). |
| GeneID | 5535 |
| Swissprot | Q96LZ3 |
| Synonyms | PPP3R2,PPP3RL |

Product Information

| | |
|----------------------|---|
| Calculated MW | 19kDa |
| Observed MW | 16kDa |
| Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purify | Affinity purification |
| Dilution | WB 1:500-1:2000 |

Background

PPP3R2 (Protein Phosphatase 3 Regulatory Subunit B, Beta) is a Protein Coding gene. Diseases associated with PPP3R2 include Arachnoiditis and Cauda Equina Syndrome. Among its related pathways are Nur77 Signaling in T-Cell and Activation of cAMP-Dependent PKA. Gene Ontology (GO) annotations related to this gene include calcium ion binding. An important paralog of this gene is PPP3R1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.