


CLPTM1 Polyclonal Antibody

Catalog No.	E-AB-62680	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using CLPTM1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human CLPTM1 (NP_001285.1).
GeneID	1209
Swissprot	O96005
Synonyms	CLPTM1

Product Information

Calculated MW	65kDa/75kDa/76kDa
Observed MW	100kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

CLPTM1 (CLPTM1 Regulator Of GABA Type A Receptor Forward Trafficking) is a Protein Coding gene. Diseases associated with CLPTM1 include Cleft Lip and Cleft Lip/Palate-Ectodermal Dysplasia Syndrome. An important paralog of this gene is CLPTM1L.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.