


UCHL5 Polyclonal Antibody

Catalog No.	E-AB-62446	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A549 cells using UCHL5 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human UCHL5 (NP_001186190.1).
GeneID	51377
Swissprot	Q9Y5K5
Synonyms	UCHL5,CGI-70,INO80R,UCH-L5,UCH37

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

UCHL5 (Ubiquitin C-Terminal Hydrolase L5) is a Protein Coding gene. Diseases associated with UCHL5 include Myeloma, Multiple. Among its related pathways are Signaling by GPCR and Metabolism of proteins. Gene Ontology (GO) annotations related to this gene include endopeptidase inhibitor activity. An important paralog of this gene is BAP1.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.