


JMJD4 Polyclonal Antibody

Catalog No.	E-AB-62516	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using JMJD4 Polyclonal Antibody at dilution of 1:5000.

Immunogen Information

Immunogen	Recombinant fusion protein of human JMJD4 (NP_075383.2).
GeneID	65094
Swissprot	Q9H9V9
Synonyms	JMJD4

Product Information

Calculated MW	50kDa/52kDa
Observed MW	60kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:1000-1:5000

Background

JMJD4 (Jumonji Domain Containing 4) is a Protein Coding gene. Diseases associated with JMJD4 include Intracranial Chondrosarcoma and Familial Isolated Hypoparathyroidism. An important paralog of this gene is JMJD6.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.