


STXBP3 Polyclonal Antibody

Catalog No.	E-AB-62935	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U2OS cells using STXBP3 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human STXBP3 (NP_009200.2).
GeneID	6814
Swissprot	O00186
Synonyms	STXBP3,MUNC18-3,MUNC18C,PSP,UNC-18C

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

STXBP3 (Syntaxin Binding Protein 3) is a Protein Coding gene. Diseases associated with STXBP3 include Hemophagocytic Lymphohistiocytosis, Familial, 4. Among its related pathways are Response to elevated platelet cytosolic Ca²⁺ and Translocation of GLUT4 to the plasma membrane. Gene Ontology (GO) annotations related to this gene include syntaxin binding and syntaxin-1 binding. An important paralog of this gene is STXBP1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.